


Aesculap Endoscopic Technology

- (GB)** Instructions for use/Technical description
- (USA)** Full HD 3CCD Camera system PV440
- (D)** Gebrauchsanweisung/Technische Beschreibung
Full HD 3CCD Kamerasystem PV440
- (F)** Mode d'emploi/Description technique
Full HD 3CCD Système de caméra PV440
- (E)** Instrucciones de manejo/Descripción técnica
Full HD 3CCD Sistema de cámara PV440
- (I)** Istruzioni per l'uso/Descrizione tecnica
Full HD 3CCD Sistema videocamera PV440
- (P)** Instruções de utilização/Descrição técnica
Full HD 3CCD Sistema da câmara PV440
- (NL)** Gebruiksaanwijzing/Technische beschrijving
Full HD 3CCD Camerasysteem PV440

- (S)** Bruksanvisning/Teknisk beskrivning
Full HD 3CCD Kamerasystem PV440
- (RUS)** Инструкция по применению/Техническое описание
Full HD 3CCD Камера PV440
- (CZ)** Návod k použití/Technický popis
Full HD 3CCD Kamerový systém PV440
- (PL)** Instrukcja użytkowania/Opis techniczny
Full HD 3CCD Układ kamery PV440
- (SK)** Návod na používanie/Technický opis
Full HD 3CCD Kamerový systém PV440
- (TR)** Kullanım Kılavuzu/Teknik açıklama
Full HD 3CCD Kamera sistemi PV440
- (KR)** 사용 설명서 / 기술 설명
의료용카메라헤드 (형명개별기재)

Aesculap®

Full HD 3CCD Camera system PV440

Legend

- 1 Endo-lens (example image)
- 2 Gray Push button
- 3 Blue Push button
- 4 Camera cable
- 5 Camera head
- 6 Camera plug
- 7 Camera control unit
- 8 ON/OFF LED
- 9 Camera display
- 10 Push button (for mode setting)
- 11 Push button (for increasing brightness)
- 12 Push button (for automatic white balance)
- 13 Camera socket
- 14 Push button (for decreasing brightness)
- 15 Push button (ON/STANDBY)
- 16 Mono jack for remote control (3.5 mm)
- 17 DVI-D video output (High Definition)
- 18 50/60 Hz switch (serial number 2000 or higher)
- 19 RGB video output (High Definition)
- 20 S-VIDEO output (Standard Definition)
- 21 VIDEO output (Standard Definition)
- 22 HD-SDI video output 1 (High Definition)
- 23 HD-SDI video output 2 (High Definition)
- 24 ON/OFF switch
- 25 Mains power
- 26 Equipotentialization connection

Symbols on product and packages

	Caution, general warning symbol Caution, see documentation supplied with the product
	Follow the instructions for use
	Marking of electric and electronic devices according to directive 2002/96/EC (WEEE). see Disposal
	Equipotentialization connection
	Fuse
	Type CF applied part
	Power on/standby
	Camera socket
	Push button for automatic white balance or to set in the camera menu
	Brightness
	Push button for mode selection or to close the camera menu
	Push button for increasing brightness or to move "up" in the camera menu
	Push button for decreasing brightness or to move "down" in the camera menu
	Low temperature and plasma sterilization

Contents

- 1. Applicable to 3
- 2. Safe handling 3
- 3. Product description 4
 - 3.1 Scope of supply 4
 - 3.2 Components required for operation 4
 - 3.3 Intended use 4
 - 3.4 Operating principle 4
- 4. Preparation and setup 4
 - 4.1 First use 4
 - Connecting the accessories 5
 - Presetting the correct voltage 5
 - Mains power 5
 - Connecting the video cable 5
 - Screwing on the endo-lens 6
 - 4.2 Settings 6
 - Access to main menu 6
 - Changing video formats 6
 - Changing the display language 6
 - Exit main menu 6
- 5. Working with the Full HD Camera 7
 - 5.1 System set-up 7
 - 5.2 Function checks 7
 - 5.3 Safe operation 7
 - Connecting the endoscope and sterile cover 7
 - Performing white balance 7
 - Selecting the mode 7
 - Adjusting light, focus and optical zoom 8
 - Adjusting the digital zoom 8
 - Adjusting the brightness 8
 - Adjusting the user mode settings 8
 - Available parameters for user modes 9
 - Remote control of an external documentation system 9
- 6. Validated reprocessing procedure 9
 - 6.1 General safety instructions 9
 - 6.2 General information 9
 - 6.3 Preparations at the place of use 10
 - 6.4 Preparation before cleaning 10
 - 6.5 Cleaning/disinfection 10
 - Product-specific safety instructions for the reprocessing procedure 10
 - Validated cleaning and disinfection procedure 10
 - 6.6 Wipe disinfection for electrical devices without sterilization 11
 - 6.7 Manual cleaning/disinfecting 11
 - Manual cleaning with immersion disinfection 11
 - 6.8 Inspection, maintenance and checks 12
 - 6.9 Packaging 12

- 6.10 Sterilization 12
- 6.11 Storage 12
- 7. Maintenance 13
- 8. Troubleshooting list 13
 - 8.1 Fuse replacement 14
- 9. Technical Service 14
- 10. Accessories/Spare parts 14
- 11. Technical data 15
 - 11.1 Ambient conditions 15
- 12. Disposal 15
- 13. Distributor in the US/Contact in Canada for product information and complaints 15

1. Applicable to

- ▶ For item-specific instructions for use and information on material compatibility, see also the Aesculap Extranet at www.extranet.bbraun.com

2. Safe handling

CAUTION

Federal law restricts this device to sale by or on order of a physician!

DANGER

Risk of fatal injury from electric shock!

- ▶ Do not open the product.
- ▶ Connect the product only to a grounded power supply.

WARNING

Risk of injury caused by incorrect operation of the product!

- ▶ Attend appropriate product training before using the product.
- ▶ For information about product training, please contact your national B. Braun/ Aesculap agency.

CAUTION

No camera image produced, due to combination of incompatible components!

- ▶ Use the camera head PV442 with serial number \geq SN3000 only with the camera control unit PV440 with serial number \geq SN2000.

- ▶ Remove the transport packaging and clean the new product, either manually or mechanically, prior to its initial sterilization.
- ▶ Prior to use, check that the product is in good working order.
- ▶ Observe "Notes on Electromagnetic Compatibility (EMC)", see TA022130.

Aesculap®

Full HD 3CCD Camera system PV440

- ▶ To prevent damage caused by improper setup or operation, and in order not to compromise warranty and manufacturer liability:
 - Only use the product according to these instructions for use.
 - Follow the safety and maintenance instructions.
 - Only combine Aesculap products with each other.
- ▶ Ensure that the product and its accessories are operated and used only by persons with the requisite training, knowledge, or experience.
- ▶ Keep the instructions for use accessible for the user.
- ▶ Always adhere to applicable standards.
- ▶ Do not operate the product in the vicinity of flammable anesthetics.
- ▶ Check the power cord for leakage currents and correct grounding at regular intervals.

3. Product description

3.1 Scope of supply

Full HD 3CCD Camera system PV440, consisting of:

- Full HD camera control unit
- Full HD 3CCD Camera head
- DVI-D cable
- Instructions for use

3.2 Components required for operation

In addition to the camera system, the following components are required for an endoscopic intervention:

- Endo-lens
- Power cord
- insert the endoscope
- Light source
- Light cable
- Monitor

3.3 Intended use

The Aesculap Full HD (=High Definition) camera is an endoscopic video camera to transmit real-time video images of a rigid or flexible endoscope, via the use of a lens with endocoupler, to a video monitor or other related documentation devices.

The Aesculap Full HD camera is suitable for all endoscopic applications.

The current use of the camera involves connecting a lens with endocoupler to the camera. An endoscope with standard eyepiece is then connected to the endocoupler and held in place by the grasping mechanism of the endocoupler.

The Full HD camera is the preferred camera for endoscopic procedures due to its superior image quality of 1920 x 1080 pixel and progressive scan mode.

3.4 Operating principle

An endoscope visualizes the operation site and an objective lens focuses this image onto the three CCD sensors, where it is converted into digital signals. The image brightness is automatically adjusted according to the light intensity (autosshutter).

Through camera cable **4**, the digital image signals are transferred to camera control unit **7**, where the digital data are processed to be displayed on the monitor.

The video signal processed in this way is then picked up at the video output monitor, and digital documentation devices.

The camera system offers automatic, electronic brightness enhancement for areas that are particularly difficult to illuminate, as well as five different fixed mode settings and three user adjustable mode settings for various situations during operation.

The video camera system is also equipped with a white adjustment system for optimal color representation.

The lens and an additional electronic zoom feature ensures that the image on the monitor is displayed at the optimum size.

4. Preparation and setup

Non-compliance with the following instructions will preclude all responsibility and liability in this respect on the part of Aesculap.

- ▶ When setting up and operating the product, adhere to
 - national regulations for installation and operation,
 - national regulations on fire and explosion protection.

Note

For the safety of patients and users it is essential that the mains power cord and, especially, the protective earth connection are intact. In many cases defective or missing protective earth connections are not registered immediately.

- ▶ Connect the device via the potential equalization terminal at the rear panel of the device to the potential equalization system of the room used for medical purposes.

Note

The potential equalization lead can be ordered from the manufacturer as art. no. GK535 (4 m length) or TA008205 (0.8 m length).

4.1 First use

WARNING

Risk of injury and/or product malfunction due to incorrect operation of the electromedical system!

- ▶ Adhere to the instructions for use of any medical device.

4.1.1 Connecting the accessories

Risk of injury due to unapproved configuration using additional components!

- ▶ For all applied components, ensure that their classification matches that of the application component (e.g. Type BF or Type CF) of the respective device.

Combinations of accessories that are not mentioned in the present instructions for use may only be employed if they are specifically intended for the respective application, and if they do not compromise the performance and safety characteristics of the products.

Also note that any equipment connected at the interfaces must demonstrably meet the respective IEC standards (e.g. IEC 60950 for data processing equipment, IEC/DIN EN 60601-1 for electromedical devices).

All configurations must comply with basic standard IEC/DIN EN 60601-1. Any individual connecting devices with one another is responsible for such configuration and must ensure compliance with basic standard IEC/DIN EN 60601-1 or applicable national standards.

- ▶ Please address your B. Braun/Aesculap partner or Aesculap Technical Service with any inquiries in this respect; for a contact address, see Technical Service.

4.1.2 Presetting the correct voltage

The mains voltage must correspond to the voltage indicated on the type plate at the back of the unit.

4.1.3 Mains power

Risk of electric shock!

- ▶ Only connect the equipment to a supply mains with protective earth.

- ▶ Plug in the power cord at mains connection **25** of camera control unit **7**.
- ▶ Connect the mains plug of the power cord to the building mains.
- ▶ Plug in the equipotentialization cable at equipotentialization connector **26** of camera control unit.
- ▶ Connect the plug of the equipotentialization cable to the building mains equipotentialization connector.
- ▶ To disconnect the camera from the mains pull out the power cord from the mains connection **25**.

4.1.4 Connecting the video cable

- ▶ Connect the video cable to the camera rear panel and to the monitor, documentation systems etc. Depending on the monitor and documentation system use the following connections:

HD Monitor (High Definition):

- DVI-D for digital HD 1080p quality
- or -
- HD-SDI for digital HD 1080i/720p quality
- or -
- RGB-Sync for analog HD 1080i/720p quality


Fig. 1

Legend

- 1 HD Monitor
- 2 SD Monitor
- 3 Full HD camera system
- 4 Mains
- 5 Equipotentialization

SD Monitor/SD documentation system (Standard Definition):

- S-Video for analog SD quality
- or -
- Video for analog SD quality


Fig. 2

Legend

- 1 HD Monitor
- 2 SD Monitor
- 3 Full HD camera system
- 4 Mains
- 5 Equipotentialization
- 6 e.g. Aesculap Eddy DVD

Aesculap®

Full HD 3CCD Camera system PV440

4.1.5 Screwing on the endo-lens

Note

Different endo-lenses are available for the camera system.


- ▶ Follow instructions for use of the endo-lens.
- ▶ Screw endo-lens 1 onto camera head 5.

4.2 Settings

4.2.1 Access to main menu

- ▶ Unplug the camera plug 6 from the camera socket 13.
- ▶ Switch on camera control unit 7.

The following test chart appears on the monitor:

- ▶ Press and hold Push button 10 for 3 seconds.

The following message appears on the camera display:

4.2.2 Changing video formats

Note

If the aspect ratio is set to 4:3, the image is cut off on the left and right margins with SD video signals.

- ▶ Select the camera frequency 50 Hz or 60 Hz (serial number 2000 and higher): Put 50/60 Hz switch 18 in the corresponding position, according to your requirements.
- ▶ Press Push button 11 in the main menu to change the video format.

The following formats can be set:

Camera setting		Video outputs		
Main menu	50/60 Hz switch	Aspect ratio 16:9		S-Video/Video
		DVI-D*	HD-SDI/RGB	
HD 1080i-SD 16:9	50 Hz	1080p 50	1080i 50	576i PAL 16:9
HD 1080i-SD 4:3				576i PAL 4:3
HD 720p-SD 16:9	50 Hz	720p 50	720p 50	576i PAL 16:9
HD 720p-SD 4:3				576i PAL 4:3
HD 1080i-SD 16:9	60 Hz	1080p 60	1080i 60	480i NTSC 16:9
HD 1080i-SD 4:3				480i NTSC 4:3
HD 720p-SD 16:9	60 Hz	720p 60	720p 60	480i NTSC 16:9
HD 720p-SD 4:3				480i NTSC 4:3

* There is always a progressive scan signal on the DVI-D output.

4.2.3 Changing the display language

- ▶ Press Push button 14 to change the display language.
 - English
 - Deutsch
 - Français
 - Español
 - Italiano

4.2.4 Exit main menu

- ▶ Press Push button 10 or plug in camera plug 6.

5. Working with the Full HD Camera

5.1 System set-up

- ▶ Position the camera tower as close as possible to the patient so that the camera cable length of 3 m can be fully utilized.
- ▶ Connect the camera plug **6** to camera socket **13** at the front panel of camera control unit **7**. Take care that the marking **UP** on the camera plug **6** is pointing upwards.

5.2 Function checks

- ▶ Check camera control unit **7** and the camera head **5** for external damage and signs of knocks or other violent impacts.
- ▶ Check that the camera cable **4** is not broken, kinked or twisted.
- ▶ Press Push button  **15** to power on the camera system.

The current settings of the camera appear on the monitor:

- ▶ Check to make certain that the front window of endo-lens **1** is completely clean. If necessary, clean the window with isopropyl alcohol (70 %).
- ▶ If necessary, ensure that the zoom adjustment ring turns freely without resistance.
- ▶ Make certain that the focusing ring turns freely without resistance.
- ▶ Align camera head **5** to an object in the room, focus the image and check for good image quality.
- ▶ To terminate the operation of the Full HD camera system press the Push button  **15**.

5.3 Safe operation

WARNING

Risk of injury and/or malfunction!

- ▶ Always carry out a function check prior to using the product.

DANGER

Risk of injury due to unapproved configuration!

- ▶ Make certain that application parts of other electrical medical devices or accessories used in this configuration for the endoscopic application are of type BF or CF.

5.3.1 Connecting the endoscope and sterile cover

- ▶ Make certain that the locking slide is opened ("Unlock" position).
- ▶ Press and hold down lock.
- ▶ Insert sterile endoscope with sterile cover and release lock.
- ▶ To lock the endoscope against rotation, close the locking slide ("Lock" position).
- ▶ Cover the camera head and cable with the sterile cover fixed to the sterile endoscope.

Note

The endocoupler can be positioned as required. To do this, swivel endocoupler to the required position!

5.3.2 Performing white balance

- ▶ Point endoscope to a white object (for example a sterile swab) at a distance of 2 cm.
- ▶ Press Push button  **12** on the control unit and keep endoscope pointed at the white object.
- or -
- ▶ Press gray Push button **2** on the camera head for longer than 2.5 seconds and keep endoscope pointed at the white object.

The following message appears on the camera display and on the monitor:

- ▶ If the white balance is performed properly, the following message appears on the camera display and on the monitor:

The image will be displayed in pure white. This procedure ensures the accurate representation of all colors.

- ▶ If the white balance did not perform properly, the following message appears on the camera display and on the monitor:

- ▶ In this case please repeat the white balance procedure.

5.3.3 Selecting the mode

To achieve the optimum image quality with endoscopes for different surgical applications, the following modes can be selected:

- Laparoscopy Standard
- Laparoscopy Advanced
- Arthroscopy
- Hysteroscopy
- Neuroendoscopy
- User 1
- User 2
- User 3

Aesculap®

Full HD 3CCD Camera system PV440

To change the mode on the camera control unit:

- ▶ Press Push button **10** on camera control unit once.
The current mode appears on the camera display and on the monitor:

- ▶ Press Push button **10** on camera control unit until the desired mode is displayed.
After 3 seconds the setting displayed will disappear and the selected mode is saved.

To change the mode on the camera head:

- ▶ Press both Push buttons on the camera head for longer than 3 seconds.
The current mode appears on the camera display and on the monitor:

- ▶ Press the blue Push button on the camera head to switch to previous mode or press the gray Push button on the camera head to switch to next mode.
After 3 seconds the setting displayed will disappear and the selected mode is saved.

5.3.4 Adjusting light, focus and optical zoom

- ▶ Adjust the brightness of the light source until sufficient illumination is achieved.
- ▶ To adjust the focus for a sharp image, turn focusing ring.
- ▶ To adjust the size of the image (zoom factor), turn the zoom adjustment ring.

5.3.5 Adjusting the digital zoom

- ▶ Press both Push buttons on the camera head and release both after a short time (2.5 seconds).
The following message appears on the camera display and on the monitor:

The following digital zoom settings can be selected:

- Zoom switched off = ZOOM 1.0
- Zoom activated = ZOOM 1.1x / 1.2x / 1.3x / 1.4x / 1.5x / 1.6x / 1.7x / 1.8x / 1.9x / 2.0x
- ▶ To change the digital zoom setting, press one of the Push buttons on the camera head as follows:
 - press the blue Push button to decrease the digital zoom.
 - press the gray Push button to increase the digital zoom. After 3 seconds the setting displayed will disappear and the selected zoom setting is saved.

5.3.6 Adjusting the brightness

To adjust the brightness on the camera control unit:

- ▶ Press Push button **11** on the camera control unit to increase the brightness.
- ▶ Press Push button **14** on the camera control unit to decrease the brightness.
The following message appears on the camera display and on the monitor:

After 3 seconds the setting displayed will disappear and the selected brightness setting is saved.

To adjust the brightness on camera head:


- ▶ Press gray Push button **2** on the camera head and release after a short time (2.5 seconds).
The following message appears on the camera display and on the monitor:

- ▶ To change the brightness, press one of the Push buttons on the camera head as follows:
 - press the blue Push button **3** to decrease the brightness.
 - press the gray Push button **2** to increase the brightness. After 3 seconds the setting displayed will disappear and the selected brightness setting is saved.

5.3.7 Adjusting the user mode settings

- ▶ Press Push button **10** on camera control unit **7** and select USER 1, USER 2 or USER 3 mode.
- ▶ Press and hold Push button **10** for 3 seconds.
The following message appears on the monitor:

- ▶ Press Push button **11** or **14** to select the desired parameter.
- ▶ Press Push button **12** to change the setting of the respective parameter.
- ▶ Press Push button **10** to save and exit user settings.
 - or –

After 40 seconds the setting options displayed will disappear and the selected settings are saved.

5.3.8 Available parameters for user modes

- Sensitivity: This parameter changes the brightness control. A low value (-) causes the camera to reduce the brightness if there is a small bright object. A high value (+) causes the camera to reduce the brightness if the average of the complete image is too bright
- Focus: This parameter decreases (-) or increases (+) the digital detail enhancement.
- Color: This parameter decreases (-) or increases (+) the color saturation.
- Gain: This parameter decreases (-) or increases (+) the maximum value of the dynamic gain. The camera reduces the dynamic gain automatically if there is enough light in the image.
- Background: This parameter decreases (-) or increases (+) the brightness of dark backgrounds.
- Red detail: This parameter decreases (OFF) or increases (LOW, HIGH) the digital enhancement especially for fine red details.
- Image size: This parameter defines the image size for the brightness measurement area. Settings are AUTO, LARGE, MEDIUM, SMALL.
- Fiberscope: The setting ON will deactivate the digital detail enhancement to minimize the moiré effect for fiberoptical endoscopes.

5.3.9 Remote control of an external documentation system

To connect an external documentation system (e.g. Aesculap Eddy DVD or any other suitable documentation system) to the mono jack for remote control **16** labeled REMOTE:

- ▶ Connect remote control cable PV968 between the camera and the documentation system.
- ▶ For remote-controlling the external documentation system, press blue Push button **3** on camera head **5**.

6. Validated reprocessing procedure

6.1 General safety instructions

Note

Adhere to national statutory regulations, national and international standards and directives, and local, clinical hygiene instructions for sterile processing.

Note

For patients with Creutzfeldt-Jakob disease (CJD), suspected CJD or possible variants of CJD, observe the relevant national regulations concerning the reprocessing of products.

Note

Mechanical reprocessing should be favored over manual cleaning as it gives better and more reliable results.

Note

Successful processing of this medical device can only be ensured if the processing method is first validated. The operator/sterile processing technician is responsible for this.

The recommended chemistry was used for validation.

Note

If there is no final sterilization, then a virucidal disinfectant must be used.

Note

For the latest information on reprocessing and material compatibility see also the Aesculap extranet at www.extranet.bbraun.com

The validated steam sterilization procedure was carried out in the Aesculap sterile container system.

6.2 General information

Dried or affixed surgical residues can make cleaning more difficult or ineffective and lead to corrosion. Therefore the time interval between application and processing should not exceed 6 h; also, neither fixating pre-cleaning temperatures >45 °C nor fixating disinfecting agents (active ingredient: aldehydes/alcohols) should be used.

Excessive measures of neutralizing agents or basic cleaners may result in a chemical attack and/or to fading and the laser marking becoming unreadable visually or by machine for stainless steel.

Residues containing chlorine or chlorides e.g. in surgical residues, medicines, saline solutions and in the service water used for cleaning, disinfection and sterilization will cause corrosion damage (pitting, stress corrosion) and result in the destruction of stainless steel products. These must be removed by rinsing thoroughly with demineralized water and then drying.

Additional drying, if necessary.

Only process chemicals that have been tested and approved (e.g. VAH or FDA approval or CE mark) and which are compatible with the product's materials according to the chemical manufacturers' recommendations may be used for processing the product. All the chemical manufacturer's application specifications must be strictly observed. Failure to do so can result in the following problems:

- Optical changes of materials, e.g. fading or discoloration of titanium or aluminum. For aluminum, the application/process solution only needs to be of pH >8 to cause visible surface changes.
- Material damage such as corrosion, cracks, fracturing, premature aging or swelling.
- ▶ Do not use metal cleaning brushes or other abrasives that would damage the product surfaces and could cause corrosion.
- ▶ Further detailed advice on hygienically safe and material-/value-preserving reprocessing can be found at www.a-k-i.org, link to Publications, Red Brochure – Proper maintenance of instruments.

Aesculap®

Full HD 3CCD Camera system PV440

6.3 Preparations at the place of use

- ▶ Remove any visible surgical residues to the extent possible with a damp, lint-free cloth.
- ▶ Transport the dry product in a sealed waste container for cleaning and disinfection within 6 hours.

6.4 Preparation before cleaning

Camera head

- ▶ Reprocess the product immediately after use.
- ▶ Use suitable cleaning/disinfecting agents if the product is put away in a wet condition.

CAUTION

Damage to the product due to improper handling when reprocessing!

- ▶ Separate the camera head from the camera control unit prior to reprocessing.
- ▶ Reprocess the camera head and the camera control unit separately.

CAUTION

Damage to, or destruction of the product caused by mechanical cleaning/disinfection!

- ▶ Only clean and disinfect the product manually.

CAUTION

Damage to the product due to inappropriate cleaning/disinfecting agents!

- ▶ Only use cleaning/disinfecting agents approved for surface cleaning. Follow the manufacturer's instructions for the respective cleaning/disinfecting agent.

CAUTION

Damage to the product due to inappropriate cleaning/disinfecting agents and/or excessive temperatures!

- ▶ Use cleaning and disinfecting agents according to the manufacturer's instructions which
 - Ensure that cleaning and disinfecting agents are approved for use on anodized aluminum, plastics (PPSU) and stainless steel.
 - do not attack softeners (e.g. in silicone).
- ▶ Observe specifications regarding concentration, temperature and exposure time.
- ▶ Do not exceed the maximum allowable temperature of 60 °C.

6.5 Cleaning/disinfection

6.5.1 Product-specific safety instructions for the reprocessing procedure

DANGER

Risk of electric shock and fire hazard!

- ▶ Unplug the device before cleaning.
- ▶ Do not use flammable or explosive cleaning or disinfecting solutions.
- ▶ Ensure that no fluids will penetrate the product.

6.5.2 Validated cleaning and disinfection procedure

Validated procedure	Special features	Reference
Wipe disinfection for electrical devices without sterilization ■ Camera control unit	–	Chapter Wipe disinfection for electrical devices without sterilization
Manual cleaning with immersion disinfection ■ Camera head	■ Suitable cleaning brush ■ 20 ml disposable syringe ■ When cleaning products with movable hinges, ensure that these are in an open position and, if applicable, move the joint while cleaning. ■ Drying phase: Use a lint-free cloth or medical compressed air (p _{max} =5 bar)	Chapter Manual cleaning/disinfecting and subsection: ■ Chapter Manual cleaning with immersion disinfection

6.6 Wipe disinfection for electrical devices without sterilization

Phase	Step	T [°C/°F]	t [min]	Conc. [%]	Water quality	Chemical
I	Wipe disinfection	RT	≥1	-	-	Meliseptol HBV wipes 50 % Propan-1-ol

RT: Room temperature

Phase I

- ▶ Remove any visible residues with a disposable disinfectant wipe.
- ▶ Wipe all surfaces of the optically clean product with a fresh, disposable disinfectant wipe.
- ▶ Observe the specified application time (1 min minimum).

6.7 Manual cleaning/disinfecting

- ▶ Prior to manual disinfecting, allow water to drip off for a sufficient length of time to prevent dilution of the disinfecting solution.
- ▶ After manual cleaning/disinfection, check visible surfaces visually for residues.
- ▶ Repeat the cleaning/disinfection process if necessary.

6.7.1 Manual cleaning with immersion disinfection

Phase	Step	T [°C/°F]	t [min]	Conc. [%]	Water quality	Chemical
I	Disinfecting cleaning	RT (cold)	>15	2	D-W	Aldehyde-free, phenol-free, and QUAT-free concentrate, pH ~ 9*
II	Intermediate rinse	RT (cold)	1	-	D-W	-
III	Disinfection	RT (cold)	15	2	D-W	Aldehyde-free, phenol-free, and QUAT-free concentrate, pH ~ 9*
IV	Final rinse	RT (cold)	1	-	FD-W	-
V	Drying	RT	-	-	-	-

D-W: Drinking water

FD-W: Fully desalinated water (demineralized, low microbiological contamination: drinking water quality at least)

RT: Room temperature

*Recommended: BBraun Stabimed

- ▶ Note the information on appropriate cleaning brushes and disposable syringes, see Validated cleaning and disinfection procedure.

Aesculap®

Full HD 3CCD Camera system PV440

Phase I

- ▶ Fully immerse the product in the cleaning/disinfectant for at least 15 min. Ensure that all accessible surfaces are moistened.
- ▶ Clean the product with a suitable cleaning brush in the solution until all discernible residues have been removed from the surface.
- ▶ If applicable, brush through non-visible surfaces with an appropriate cleaning brush for at least 1 min.
- ▶ Mobilize non-rigid components, such as set screws, links, etc. during cleaning.
- ▶ Thoroughly rinse through these components with the cleaning disinfectant solution (at least five times), using a disposable syringe.

Phase II

- ▶ Rinse/flush the product thoroughly (all accessible surfaces) under running water.
- ▶ Mobilize non-rigid components, such as set screws, joints, etc. during rinsing.
- ▶ Drain any remaining water fully.

Phase III

- ▶ Fully immerse the product in the disinfectant solution.
- ▶ Mobilize non-rigid components, such as set screws, joints, etc. during rinsing.
- ▶ Rinse lumens at least 5 times at the beginning of the exposure time using an appropriate disposable syringe. Ensure that all accessible surfaces are moistened.

Phase IV

- ▶ Rinse/flush the product thoroughly (all accessible surfaces).
- ▶ Mobilize non-rigid components, such as set screws, joints, etc. during final rinse.
- ▶ Rinse lumens with an appropriate disposable syringe at least five times.
- ▶ Drain any remaining water fully.

Phase V

- ▶ Dry the product in the drying phase with suitable equipment (e.g. cloth, compressed air), see Validated cleaning and disinfection procedure.

6.8 Inspection, maintenance and checks

- ▶ Allow the product to cool down to room temperature.
- ▶ Inspect the product after each cleaning and disinfecting cycle to be sure it is: clean, functional, and undamaged.
- ▶ Check the product for any damage, abnormal running noise, overheating or excessive vibration.
- ▶ Set aside the product if it is damaged.
- ▶ Remove the moisture from the camera plug and check that the contacts are dry.

6.9 Packaging

- ▶ Follow the instructions for use for the applied packaging and storage systems (e.g. instructions for use TA009721 for Aesculap Eccos storage system).
- ▶ Insert the product in its proper position in the Eccos holder, or put it on a tray in such a way that the product is protected against damage. Ensure that all cutting edges are protected.
- ▶ Pack trays appropriately for the sterilization process (e.g. in Aesculap sterile containers).
- ▶ Ensure that the packaging will prevent a recontamination of the product.

6.10 Sterilization

CAUTION

Damage to the camera control unit by sterilization!

- ▶ Do not sterilize the camera control unit under any circumstances.
- ▶ Clean the camera control unit by surface disinfection only.

CAUTION

Damage to the camera head or the optical zoom coupler caused by inappropriate sterilization processes!

- ▶ Never sterilize the camera head or the optical zoom coupler with steam or in an autoclave.
- ▶ Do not carry out chemical sterilization.

Working with disposable sterile covers

- ▶ Use suitable Aesculap disposable sterile covers. Follow the instructions for use of the sterile covers.

Sterilization with Sterrad® 50, 100, 100S, 200 sterilization process

Note

The Sterrad® sterilization process can cause cosmetic changes to the camera head and cable. These changes, however, do not affect its functionality.

- ▶ Sterilize by the Sterrad® sterilization process 50, 100, 100S or 200, observing the following rules: Follow the manufacturer's instructions for use of the Sterrad® system. Use of a biological indicator is recommended to confirm effective sterilization.

6.11 Storage

- ▶ Store sterile products in germ-proof packaging, protected from dust, in a dry, dark, temperature-controlled area.

7. Maintenance

To ensure reliable operation, the product must be maintained at least once a year.

For technical service, please contact your national B. Braun/Aesculap agency, see Technical Service.

8. Troubleshooting list

Malfunction	Cause	Remedy
No image on the monitor	Device not powered	Verify that all units are connected to mains connector and switched on
	Fuse blown	Replace camera fuse T0.5 A L/250 V~ (TA020394)
	Wrong video input activated	If the monitor supports several channels: Make certain the monitor is switched to the channel that is connected to the camera system
	Video cable defective	Replace video cable
	Video cable not connected	Connect video cable
Image too dark	Camera cable not connected or defective	Ensure that the camera head is connected to the camera control unit and the camera plug is dry
	Light cable not connected	Verify that the light cable is connected at the light source and the light cable input of the endoscope
	Light cable defective	Connect a new light cable/Have repaired by manufacturer
	Light source dimmed down	Turn up the light source intensity
	Endoscope optical system defective	If the image is too dark even without the camera: Use a different endoscope and have the defective endoscope repaired by the manufacturer
Image fuzzy or poorly defined	Object out of focus	Turn the focus adjustment setting ring of the zoom coupler until there is a sharp image
	Debris on the endoscope tip	Clean tip with a sterile cloth and isopropyl alcohol (70 %)
	Moisture on the camera cable plug	Dry plug with gauze before plugging it in at the camera plug
Image fogged	Contaminated or fogged endoscope	Check if the front and rear windows of the endoscope and the front window of the zoom coupler are clean; if necessary, clean with isopropyl alcohol (70 %) Use Aesculap Anti Fog Solution JG910 on the window on the distal end of the endoscope
Noise or excessive grain in the picture	Brightness enhancement function is active	Increase intensity of the light source
Image too bright or glary	Light metering not at its optimum setting	Press Push button 10 for the mode selection and select the optimum light metering mode

Aesculap®

Full HD 3CCD Camera system PV440

Malfunction	Cause	Remedy
Color bar displayed on the monitor	Camera cable not connected to the camera control unit	Connect camera cable at the camera control unit
Inaccurate color representation on the monitor	White balance failed	Initiate white-adjustment
	Inappropriate monitor color settings	Reset to factory setting
	Video cable defective	Replace video cable

8.1 Fuse replacement

Risk of fatal injury from electric shock!
 ► Unplug the device before changing the fuses!

Specified fuses: T 0.5 A L/250 V~; Art. no. TA020394

- Use a small screwdriver to release the clip on the fuse holder.
- Pull out the fuse holder.
- Replace both fuses.
- Reinsert the fuse holder in such a way that it audibly snaps into place.

Fig. 3 Fuse replacement

Note

If the fuses burn out frequently, the device is faulty and should be repaired, see Technical Service.

9. Technical Service

Risk of injury and/or malfunction!
 ► Do not modify the product.

- For service and repairs, please contact your national B. Braun/ Aesculapagency.

Modifications carried out on medical technical equipment may result in loss of guarantee/warranty rights and forfeiture of applicable licenses.

Service addresses

Repairs handling camera systems:

Aesculap FLEXIMED GmbH

Repair Service

Robert-Bosch-Str. 1

79211 Denzlingen/Germany

Phone: +49 7666 93 21 0

Fax: +49 7666 93 21 580

E-mail: techserv@aesculap-fleximed.de

Telephone hotline for technical assistance with electrical devices and device combinations:

Phone: +49 7461 95 1601

Other service addresses can be obtained from the address indicated above.

10. Accessories/Spare parts

Art. no.	Designation
PV442	Camera head and camera cable
PV441	Camera cable
PV126S	Endo-lens with zoom
PV127S	Fixed focal length endo-lens
TA020394	Fuse set T 0.5 A L/250 V~

11. Technical data

Classification acc. to Directive 93/42/EEC

Art. no.	Designation	Class
PV442	Camera head and camera cable	I
PV440	Camera control unit	

Art. no.	Designation
Mains power	100-240 V AC, 50 Hz to 60 Hz
Power consumption	45 VA
Image sensors	3 CCD sensor 1/3 Inch
Resolution	1920 x 1080 Pixel
Electronic Zoom	1.0x to 2.0x
Video outputs:	
DVI-D	1080p or 720p*; 50 Hz or 60 Hz** Aspect ratio 16:9
HD-SDI	1080i or 720p*; 50 Hz or 60 Hz** Aspect ratio 16:9
RGB	1080i or 720p*; 50 Hz or 60 Hz** Aspect ratio 16:9
S-Video	PAL or NTSC** Aspect ratio 16:9 or 4:3*
Video	PAL or NTSC** Aspect ratio 16:9 or 4:3*
Protection class (acc. to IEC/DIN EN 60601-1)	Camera head IPX7 Control unit IP20
Dimensions of control unit (L x W x H)	305 mm x 75 mm x 305 mm
Dimensions of camera head (without zoom lens) (L x W x H)	40 mm x 53 mm x 72 mm
Applied part	Type CF
EMC	IEC/DIN EN 60601-1-2
Conforming to stan- dard	IEC/DIN EN 60601-1

* Video format configurable in the main menu, see Changing video formats

** 50 Hz/60 Hz or PAL/NTSC are switchable on the rear panel of the device from series number 2000 and higher

11.1 Ambient conditions

	Operation	Storage and transport
Temperature	+10 °C to +40 °C	-20 °C to +60 °C
Relative humidity	0 % to 90 % (at +20 °C) or 0 % to 50 % (at +40 °C), non-condensing	0 % to 90 % (at +20 °C) or 0 % to 50 % (at +40 °C), non-condensing

12. Disposal

Note

The user institution is obliged to process the product before its disposal, see Validated reprocessing procedure.

	<p>Adhere to national regulations when disposing of or recycling the product, its components and its packaging! The recycling pass can be downloaded from the Extranet as a PDF document under the respective article number. (The recycling pass includes disassembling instructions for the product, as well as information for proper disposal of components harmful to the environment.) Products carrying this symbol are subject to separate collection of electrical and electronic devices. Within the European Union, disposal is taken care of by the manufacturer as a free-of-charge service.</p>
---	---

- Detailed information concerning the disposal of the product is available through your national B. Braun/Aesculap agency, see Technical Service.

13. Distributor in the US/Contact in Canada for product information and complaints

3773 Corporate Parkway
Center Valley, PA, 18034,
USA

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Aesculap AG | Am Aesculap-Platz | 78532 Tuttlingen | Germany