STORZ

KARL STORZ — ENDOSKOPE

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20 113320 Kaltlicht-Fontäne HALOGEN 250 twin

INSTRUCTION MANUAL

20 113320 Cold light fountain HALOGEN 250 twin

MANUAL DE INSTRUCCIONES

20 113320 Fuente de luz fría HALOGEN 250 twin

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**Important information
for users of KARL STORZ
instruments**

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It is recommended to check the suitability of the product for the intended procedure prior to use.

Thank you for your expression of confidence in the KARL STORZ brand name. Like all of our other products, this product is the result of years of experience and great care in manufacture. You and your organization have decided in favor of a modern, high-quality item of equipment from KARL STORZ.

This instruction manual is intended to serve as an aid in the proper setup, installation, and operation of the HALOGEN 250 twin cold light fountain. All essential details of the equipment and all actions required on your part are clearly presented and explained. We thus ask that you read this manual carefully before proceeding to work with the equipment. Insert this manual in its protective wallet and keep it available for ready reference in a convenient and conspicuous location near the equipment.

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


	Read the instructions carefully before operating the equipment!
	Lamp defective indicator
	Pilot lamp
	Luminous intensity switch
	Light exit
	Applied part of type CF
	Potential equalization connector
	Lamp fuse
	Warning: Burn hazard! Lamp and its surroundings are hot.
	Lamp type (Art. no. 105)
	Electronic information product pollution control (China RoHS)
	Manufacturer

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 **DANGER:** Risk of explosion if used in the presence of flammable anesthetics.

 **CAUTION:** To reduce the risk of electrical shock, do not remove cover. Refer servicing to qualified service personnel.

Keep out of reach of patients.
Do not store liquids on or above the unit.

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Device description

The cold light fountain HALOGEN 250 twin is suitable for all types of endoscopic examination. Color temperature is approx. 3.400 K.

A high-quality condenser system ensures high luminous efficiency while modern design guarantees easy operation and unproblematic servicing.

In case of possible failure of the main lamp, a second halide lamp, kept in reserve, is automatically put into operation. This guarantees that, if a lamp does fail, the operation may be continued at the same high lighting power without interruption.

Adaptors for light cables of other manufacturers are available as accessories.

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Warnings and cautions

Please read this manual and follow its instructions carefully. The words **Warning**, **Caution**, and **Note** convey special meanings. Wherever they are used in this manual, they should be carefully reviewed to ensure the safe and effective operation of this product. To make the words stand out more clearly, they are accompanied by a pictogram.

WARNING: A Warning indicates that the personal safety of the patient or physician may be involved. Disregarding a Warning could result in injury to the patient or physician.

CAUTION: A Caution indicates that particular service procedures or precautions must be followed to avoid possible damage to the product.

NOTE: A Note indicates special information about operating the product, or clarifies important information.

WARNING: Read this instruction manual thoroughly and be familiar with its contents prior to using this equipment. Before using the unit, read the following safety instructions carefully to avoid putting your patients, personnel or yourself at risk.

NOTE: Any damage to the unit resulting from incorrect operation is not covered by the manufacturer's warranty.

NOTE: Lamps do not have an indefinite service life, so it is possible that one may fail during use. For this reason a spare lamp system or cold light source should always be kept available during therapeutic endoscopic interventions.

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WARNING: The electrical installations in the operating room in which the unit is connected and operated must comply with the applicable IEC standards.

WARNING: Keep out of reach of patients.

WARNING: The instructions and interface specifications for medical devices and/or system components used in combination must be observed precisely.

WARNING: Grounding reliability can only be achieved when the equipment is connected to a properly installed "Hospital Only" or "Hospital Grade" receptacle (i.e., approved for use in an operating room environment). Routinely inspect electrical plug and cord. Do not use if inspection reveals damage.

WARNING: Test this equipment prior to each surgical procedure to ensure that it functions correctly.

WARNING: Risk of explosion if used in the presence of flammable anesthetics.

WARNING: Always unplug the unit before doing any maintenance work on it (e.g. cleaning).

WARNING: Danger of burns. Lamp may be hot.

WARNING: Light cable socket may be very hot; particularly if the temperature switch has been triggered and the unit has switched off.

WARNING: Danger of electric shock. Do not open the unit. Refer servicing only to qualified personnel. Any opening of the unit by unauthorized persons voids the guarantee. The only exception to this are activities described in this instruction manual.

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CAUTION: Always adjust the cold light fountain to the minimum illumination intensity necessary to achieve optimum illumination of the endoscopic scene, either by direct vision or when coupled to a video camera.

CAUTION: Avoid allowing fluids to enter the unit. Do not store liquids on or above the unit.

CAUTION: Before connecting the unit to the electrical supply lines, verify that the line voltage stated on the unit's identification plate agrees with that of the electrical supply lines to be used.

CAUTION: Use only fuses of the correct rating.

CAUTION: Do not touch glass bulb with your bare fingers, use the protective cover. Touching the lamp with bare fingers may result in starting difficulty.

CAUTION: Federal law restricts this device to sale by or on the order of a physician (USA only).

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Light cables

WARNING: Danger of glare! Never look into the open end of a light cable.

WARNING: Danger of glare! Do not operate the light source without a light cable adapter.

WARNING: Danger of burns. Contact with the open end of a light cable connected to a cold light fountain can possibly lead to burns.

WARNING: High energy radiated light through endoscopes may give rise to high temperatures in front of the light outlet and at the tip of the endoscope.

WARNING: Never leave an illuminated light cable with or without an endoscope connected, on or near a patient; the risk of patient burns can result.

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CAUTION: Do not place connected light cables on combustible materials such as textiles (drapes) or in the direct vicinity of swabs, etc. soaked with combustible fluids, such as disinfectant. These could catch fire from heat development at the open end of the light cable. Therefore always attach the open end of a connected light cable to the provided retainer ⑤.

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Before using the unit on the patient it is imperative that you be acquainted with how the unit operates and is controlled.

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Intended use

The cold light fountain HALOGEN 250 twin model 20 113320 is used for illumination in diagnostic and therapeutic applications in endoscopy.

Use of the unit in fields other than those indicated above is not allowed for safety reasons.

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The cold light fountain HALOGEN 250 twin may only be used with accessories which have been designated by KARL STORZ as suitable for the unit or the safe use of which is proven.

i NOTE: KARL STORZ light cables are recommended for use with this unit. Some manufacturer's cables may be damaged when this unit is operated at higher settings.

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Unauthorized conversions or modifications to the unit are not allowed for safety reasons.

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
User qualification

The cold light fountain HALOGEN 250 twin may only be used by physicians and medical assistants who have a corresponding specialized qualification and who have been instructed in use of the unit.

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Safety precautions at the site of installation

The unit may only be used in medical rooms in which the electrical equipment has been installed in accordance with applicable national regulations.

It is not intended for use in hazardous zones. This means, for example, that when using easily combustible and explosive inhalation anesthetics or mixtures thereof, the unit must not be operated inside the hazard zone shown in the diagram. This also applies to highly combustible and potentially explosive chemicals, e.g. skin disinfectants and fast-acting surface disinfectants.

The unit is equipped with a connector for attaching a ground line. It should be connected up in accordance with national regulations.

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Safety precautions when operating the unit

It is the user's responsibility to make sure the unit is safe and operates properly before using the unit.

During treatment using the cold light fountain HALOGEN 250 twin the patient must be treated and kept under observation with the usual medical care. This includes keeping a check on the progress of treatment, as well as monitoring sterile application conditions where required by the type of intervention.

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Unpacking the equipment

Carefully unpack the cold light fountain and remove the unit and its accessories from their packaging.

Check for missing items and evidence of shipping damage.

File any complaints with the manufacturer or supplier immediately. If possible, retain the original packing materials for later use; these can come in handy if the unit has to be transported.

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
Basic equipment

- 1 Cold light fountain HALOGEN 250 twin
- 1 Power cord 400 A
- 1 Instruction manual
- 1 Wallet for instruction manual

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Installation

! **NOTE:** The cold light fountain HALOGEN 250 twin may only be used in medical rooms in which the electrical equipment has been installed in accordance with applicable national regulations.

! **WARNING:** This unit is not intended for use in hazardous zones. Do not operate the unit within the hazard zone shown in the diagram while explosive anesthetic gases are in use.

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Set the unit on a flat surface.

! **CAUTION:** Do not cover the air vent. Danger of overheating. Make certain there is free air circulation towards the top and to the side (distance of rear panel to wall: 10 cm). If the unit overheats it will shut down automatically. Turn off the main power switch ① and let the light source cool down (see error list).

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The unit is equipped with a connector ⑦ for attaching a ground line.

The unit's ground line should be installed by a qualified electrician.

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
**Operating the
cold light fountain**
Setting of the system voltage

Remove the power fuse holder ⑨ (see section "Fuse replacement", page 16) and take out the voltage selection module.

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Insert the voltage selection module into the power fuse holder in such a way that the system voltage displayed in the viewing window on the back of the power fuse holder corresponds to the local power supply.

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
Check the inserted power fuses using the following table:

100 VAC	2 x T6.3 AL250V
120 VAC	2 x T6.3 AL250V
230 VAC	2 x T3.15 AL250V
240 VAC	2 x T3.15 AL250V

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Connect power cord. Insert power cord into power cord receptacle ⑩ as far as it will go.

! WARNING: Only insert the power plug into and remove it from electrical outlets located outside areas subject to explosion hazards.

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Plug the light cable into socket ⑥ as far as it will go. To disconnect the light cable pull on the connector, never on the cable itself.

! WARNING: Danger of glare. Never look into the open end of a light cable.

! WARNING: Danger of burns. Contact with the open end of a light cable connected to a cold light fountain can possibly lead to burns.

! WARNING: Light cable socket may be very hot; particularly if the temperature switch has been triggered and the unit has switched off.

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Connect the light cable to the endoscope by a quarter turn of the knurled collar on the screw base.

i NOTE: KARL STORZ light cables are recommended for use with this unit. Some manufacturer's cables may be damaged when this unit is operated at higher settings.

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If you do not plan to conduct the light cable to the endoscope right away, plug the loose end into the dummy socket ⑤.

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① NOTE: The entrance and exit faces of the prepared light cable and endoscope must be kept clean if optimal light delivery is to be maintained. Remove crude stains with a lint-free disposable cloth, disinfect the entrance and exit faces with a cotton tip applicator soaked in alcohol and dry.

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Switch on the unit (operating pilot lamp ③ is lit).

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① NOTE: The light exit ⑥ is provided with an antiglare flap which does not allow any direct outlet of light. Light only appears when the light cable is attached.

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Adjust the required luminous intensity with the intensity switch ④ (3 levels).

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⚠ WARNING: Danger of glare. Never look into the open end of a light cable.

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⚠ WARNING: Danger of burns. Contact with the open end of a light cable connected to a cold light fountain can possibly lead to burns.

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i NOTE: Frequently switching the cold light fountain on and off causes the halogen lamp to wear out more quickly. For relatively brief interruptions, set the luminous intensity switch ④ to Level 1.

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If the halogen lamp is defective, the lamp defective indicator ② lights red and, at the same time, the green operating indicator ③ goes out. The device switches to the spare lamp. The defective lamp should be replaced as soon as possible (see section 'Lamp replacement', page 13).

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Adaptor for light cables of other manufacturers

If necessary light exit ⑥ can be unscrewed and replaced by adaptors for light cables of other manufacturers (see KARL STORZ catalog TELEPRESENCE).

To do so, unscrew light exit ⑥ by rotating counterclockwise.

! WARNING: Danger of glare! Do not operate the light source without a light cable adapter.

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Place the adaptor on the screw thread and tighten up by rotating clockwise.

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
Lamp replacement

Switch off the unit. Disconnect the power plug.

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Open the housing cover. For this purpose, first undo the 2 thumbscrews at the bottom of the unit.

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Turn the unit back up the right way and slide the housing cover towards the back.

i NOTE: Slide back the housing cover just far enough that the lamp compartment is easily accessible.

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
Remove the lid from the lamp compartment. To do this, loosen the thumbscrew.

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i NOTE: If there is a lamp defect, the unit switches over automatically to the spare lamp. To do this, the position of the lamp holder is changed mechanically (The spare lamp swings into the main lamp position). After replacing the lamp, the lamp base must be returned counterclockwise to the "normal position" (until it clicks into place) using the accompanying hex key. The main lamp is directly in front of the reflector.

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Pull out the defective halogen lamp (left plug-in position, see arrow) and discard.

! WARNING: Danger of burns. Lamp may be hot.

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
Remove the spare lamp from the right plug-in position.

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Replace the main lamp with the spare lamp

i NOTE: The plug-in position for the main lamp is more easily accessible if the lamp base is in the normal position.

Using the hex key, turn the lamp base counterclockwise into the normal position "until it clicks into place" (see above sketch).

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Insert the spare lamp into the plug-in position for the main lamp (left) (use the protective cover!).

! CAUTION: When replacing the lamps, check for easy movement of the lamp replacement unit (lamp base). There should be nothing in the lamp compartment which could block the turning motion of the lamp base.

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Inserting a new lamp in the plug-in position for the spare lamp (right)

! **CAUTION:** Do not touch the glass of the new bulb with your fingers.

i **NOTE:** The plug-in position for the spare lamp is more easily accessible if the lamp base is in the "defect position."


If necessary, turn the lamp base to the "defect position" using the hex key.

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i **NOTE:** Do not touch glass bulb with your bare fingers. (If necessary, wipe the bulb with a cotton swab moistened with alcohol before operation.)

Please use the protective cover provided with the new lamp. Insert new halogen lamp (24 V/250 W) into the contact socket for the plug-in position of the spare lamp until firmly in place (see arrow). The lamp is automatically adjusted.

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i **NOTE:** The former spare lamp is now used as the main lamp, while the spare lamp is replaced by a new halogen lamp. This ensures that a "new" spare lamp is always available.

! **CAUTION:** Before assembling the cold light fountain, make certain that the lamp base has been returned to the normal position.

If necessary, turn the lamp base to the normal position.

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Reinstall the lamp compartment lid.
Refit the housing cover and tighten the thumbscrews.
Reconnect the power cord.
Carry out function test.

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Fuse replacement

Switch off the unit and remove the power cord.

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
Remove the power fuse holder ⑨ with a suitable tool.

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! **CAUTION:** Only use fuses of the correct rating.

Insert new fuses of the appropriate rating.

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Refit the power fuse holder ⑨.
Reconnect the power cord.
Test the instrument for proper operation.

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Preparation

WARNING: Always pull out power plug before cleaning! Make sure the lamp has cooled down.

CAUTION: Avoid allowing fluids to enter the unit.

The exterior surfaces of the medical device must be wiped clean with a disposable cloth moistened with disinfectant or a ready to use soaked disinfectant cloth. Due to their protein-fixating effect and material incompatibility, alcohol-based agents must not be used.

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
Maintenance and safety check

Maintenance

Preventive maintenance is not essential. Regular maintenance can, however, contribute to identifying potential problems before they become serious, thus enhancing the instrument's reliability and extending its useful service life. Maintenance services can be obtained from your local representative or from the manufacturer.

Safety checks/repeat inspections as per IEC 62353

WARNING: Regardless of the accident prevention regulations or testing intervals for medical devices prescribed in different countries, with this device safety checks/repeat inspections as defined by IEC 62353 must be performed and recorded by a qualified electrician at least once a year.

Visual inspection

1. Check the device and accessories for any mechanical damage which may impair functionality.
2. Check that inscriptions relevant to safety are legible.

Electric measurements

- Inspect the device safety fuses
- Measure protective ground resistance in accordance with IEC 62353: The limit values can be found in the current standard.
- Measure earth leakage current in accordance with IEC 62353: The limit values can be found in the current standard.
- Touch current measured according to IEC 62353: The limit values can be found in the current standard.
- Patient leakage current measured according to IEC 62353: The limit values can be found in the current standard.

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Test for proper functioning

Carry out the functional check as per the instruction manual (see section "Operation").

Documentation

The safety check/repeat inspection and results must be documented.

WARNING: If any defects or shortcomings which could endanger patients, personnel or others are ascertained during these safety checks/repeat inspections, the device must not be operated until such time as these defects or shortcomings have been eliminated by qualified technical servicing.

NOTE: For detailed information on the range and implementation of the safety check/repeat inspection, please refer to the latest issue of the relevant Service Manual.

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Servicing and repair

Defective items of equipment must be serviced and repaired exclusively by persons authorized by KARL STORZ; all repair work must employ original KARL STORZ parts only.

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Disposal

This unit has been marked in accordance with the European Directive on waste electrical and electronic equipment (WEEE).

At the end of its useful operating life, dispose of the unit as electronic scrap.

Please ask either KARL STORZ SE & Co. KG, a KARL STORZ subsidiary or your specialist dealer for information on your local collection point.

Within the scope of application of this Directive, KARL STORZ SE & Co. KG is responsible for the proper disposal of this unit.

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Repair program

Usually to bridge the repair period, you will receive a unit on loan which you then return to KARL STORZ as soon as you receive the repaired unit.

In Germany you can refer repairs directly to

KARL STORZ SE & Co. KG
Abt. Reparaturservice
Take-off Gewerbepark 83
78579 NEUHAUSEN

In other countries please contact your local KARL STORZ branch or authorized dealer.

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Important information

In order to prevent infection, it is strictly forbidden to ship contaminated medical devices. All medical devices must be decontaminated on site, to avoid contact and aerogenous infections (among staff).

We reserve the right to return contaminated instruments/equipment to the sender.

Repairs, changes or expansions which are not carried out by KARL STORZ or by experts authorized by KARL STORZ will invalidate all guarantee rights. KARL STORZ gives no guarantee on the correct functioning of equipment or instruments which have been repaired by unauthorized third parties.

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Limitation of liability

KARL STORZ shall be liable for the safe operation, operational reliability, and performance of this equipment only subject to the following conditions:

- all assembly operations, system expansions, readjustments, modifications, or repairs have been performed by a person or persons duly authorized by KARL STORZ;
- all electrical installations at the location in which the unit is connected and operated meet the applicable laws and standards; and
- the unit has been used in accordance with its operating instructions at all times.

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Manufacturer's warranty

The guarantees provided can be found in the Standard Conditions of Business of KARL STORZ.

The medical device must always be sent to your local subsidiary (see "Subsidiaries" page), even during the warranty period.

Opening the equipment or performance of any repairs or modifications to the equipment by unauthorized persons shall relieve us of any liability for its performance. Any such repair or modification performed during the warranty period shall void all warranty.

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Troubleshooting

WARNING: Always unplug the unit before all maintenance work.

Symptom:
Complete failure of the unit.

Possible causes:

- No power from the power line.
- Defective power fuse.
- Temperature safety device triggered.

Remedy:

- Check that there is electricity to the wall outlet.
- Change fuses as described in the instruction manual. Make sure you use the correct type of fuse.
- Disconnect unit from power and allow to cool down. Switch on again after approx. 10–15 minutes. Make certain there is sufficient ventilation.

Symptom:
No light is emitted, red pilot lamp lights.

Possible causes:
Spare lamp is defective.

Remedy:
Insert new lamp.

Symptom:
Too little light.

Possible causes:

- Luminous intensity switch is set too low (Stage 1 or 2).
- Ends of light cable and/or endoscope are dirty.
- Light cable is defective.
- Dimmer defective.

Remedy:

- Set the luminous intensity switch to a higher stage (2 or 3)
- Clean ends of light cables, and light entrances and exits on the endoscope.
- Replace light cable.
- Have the unit repaired.

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Symptom:
No light output in spite of changing the bulb.

Possible causes:

- Bulb inserted is defective.
- Automatic cutout (I) has operated.
- Bulb base is defective.

Remedy:

- Fit a new bulb.
- Activate the automatic cutout (I) by pressing.
- Have the unit repaired.

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Standard compliance (for 20 1133 20)

According to IEC 60601-1:

- Type of protection against electric shocks: Protection Class I
- Degree of protection against electric shocks: Applied part of type CF

According to IEC 60601-1-2:

Please read the Electromagnetic Compatibility Information in the appendix (p. 28-39).

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Directive compliance (for 20 1133 20)

According to Medical Device Directive (MDD):

This medical product belongs to Class I

This medical product bears the CE mark in accordance with MDD 93/42/EEC.

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Technical documentation

On request, the manufacturer will provide those circuit diagrams, itemized parts listings, descriptions, sets of adjustment instructions, and other items of available documentation to suitably qualified user personnel duly authorized by the manufacturer for their use in repairing those components of the instrument that have been designated by their respective manufacturers as repairable.

Supply of such technical documentation relating to the instrument shall not be construed as constituting manufacturer's authorization of user's personnel, regardless of their level of technical training, to open or repair the instrument.

Explicitly exempted herefrom are those maintenance and repair operations described in this manual.

We reserve the right to make engineering modifications in the interest of promoting technological progress and generating performance improvements without obligation on our part to submit prior notice thereof.

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Appendix

Electromagnetic Compatibility (EMC) Information

WARNING: Medical electrical equipment needs special precautions regarding Electromagnetic Compatibility (EMC). Observe the EMC instructions in this appendix during installation and commissioning.

The HALOGEN 250 twin model 20 113320 corresponds to IEC 60601-1-2 [CISPR 11 Class B] and therefore meets the EMC requirements of the Medical Device Directive (MDD) 93/42/EEC.

These limits are designed to provide reasonable protection against the typical electromagnetic interference to be expected in a medical environment. The HALOGEN 250 twin model 20 113320 is a Group 1 unit (as per CISPR 11). Group 1 contains all the 'equipment and systems which generate or use RF energy only for their internal functioning'.

NOTE: The tables and guidelines that are included in this appendix provide information to the customer or user that is essential in determining the suitability of the equipment or system for the electromagnetic environment of use, and in managing the electromagnetic environment of use to permit the equipment or system to perform its intended use without disturbing other equipment and systems or non-medical electrical equipment. If this equipment does cause harmful interference with other devices, the user is encouraged to try to correct the interference by one or more of the following measures:

- reorient or relocate the receiving device
- increase the separation between the equipment
- Connect the equipment into an outlet on a circuit different from that to which the other device(s) is connected.

If you have any further questions, please contact your local representative or our service department.

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WARNING: The HALOGEN 250 twin model 20 1133 20 should not be used adjacent to or stacked with other equipment. If adjacent or stacked use is necessary, the equipment or system should be observed to verify normal operation in the configuration in which it will be used.

WARNING: The use of portable and mobile RF equipment may have an impact on this and other pieces of medical equipment.

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Table 201 Guidance and manufacturer's declaration – electromagnetic emissions		
<p>HALOGEN 250 twin model 20 1133 20 is intended for use in the electromagnetic environment specified below. The customer or user of the HALOGEN 250 twin model 20 1133 20 should ensure that it is used in such an environment.</p>		
Emissions test	Compliance	Electromagnetic environment – guidance
RF emissions CISPR 11	Group 1	HALOGEN 250 twin model 20 1133 20 uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause interference in nearby electronic equipment.
RF emissions CISPR 11	Class B	The HALOGEN 250 twin model 20 1133 20 is suitable for use in all establishments (e.g. hospitals and doctors' offices) including domestic establishments and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes.
Harmonic emissions IEC 61000-3-2	Class A	
Voltage fluctuations/flicker emissions IEC 61000-3-3	Complies	

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Table 202			
Guidance and manufacturer's declaration – electromagnetic immunity			
<p>HALOGEN 250 twin model 20 1133 20 is intended for use in the electromagnetic environment specified below. The customer or user of the HALOGEN 250 twin model 20 1133 20 should ensure that it is used in such an environment.</p>			
Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment – guidance
Electrostatic Discharge (ESD) IEC 61000-4-2	± 6 kV contact ± 8 kV air	Complies ± 6 kV contact ± 8 kV air	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30 %.
Electrical fast transient/burst IEC 61000-4-4	± 2 kV for power supply lines ± 1 kV for input/output lines	Complies ± 2 kV for power supply lines ± 1 kV for input/output lines	Mains power quality should be that of a typical commercial or hospital environment.
Surge IEC 61000-4-5	± 1 kV differential mode ± 2 kV common mode	Complies ± 1 kV differential mode ± 2 kV common mode	Mains power quality should be that of a typical commercial or hospital environment.
Voltage dips, short interruptions and voltage variations on power supply input lines IEC 61000-4-11	<5 % U_T * (>95 % dip in U_T) for 0.5 cycle 40 % U_T (60 % dip in U_T) for 5 cycles 70 % U_T (30 % dip in U_T) for 25 cycles <5 % U_T (>95 % dip in U_T) for 5 sec	Complies <5 % U_T * (>95 % dip in U_T) for 0.5 cycle 40 % U_T (60 % dip in U_T) for 0.5 cycle See guidelines 70 % U_T (30 % dip in U_T) for 1 cycle See guidelines Complies <5 % U_T (>95 % dip in U_T) for 5 sec	Mains power quality should be that of a typical commercial or hospital environment. If the user of the equipment or system requires continued operation during power cuts, short term power interruptions or if power supply is erratic, it is recommended that the equipment or system be powered from an uninterruptible power supply.
Power frequency (50/60 Hz) magnetic field IEC 61000-4-8	3 A/m	Complies 3 A/m	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.
* Note: U_T is the a.c. mains voltage prior to application of the test level.			

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Table 204			
Guidance and manufacturer's declaration – electromagnetic immunity – for equipment and systems that are not life-supporting			
<p>HALOGEN 250 twin model 20 1133 20 is intended for use in the electromagnetic environment specified below. The customer or user of the HALOGEN 250 twin model 20 1133 20 should ensure that it is used in such an environment.</p>			
Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment – guidance
<p>Conducted RF IEC 61000-4-6</p> <p>Radiated RF IEC 61000-4-3</p>	<p style="text-align: center;">$3 V_{rms}$ 150 kHz to < 80 MHz</p> <p style="text-align: center;">$3 V/m$ 80 MHz to 2.5 GHz</p>	<p style="text-align: center;">$3 V_{rms}$</p> <p style="text-align: center;">$3 V/m$</p>	<p>Portable and mobile RF communications equipment should be used no closer to any part of the HALOGEN 250 twin model 20 1133 20, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter.</p> <p>Recommended separation distance:</p> <p>$d = [3.5/3]\sqrt{P}$ 150 kHz to < 80 MHz</p> <p>$d = [3.5/3]\sqrt{P}$ 80 MHz to < 800 MHz</p> <p>$d = [7/3]\sqrt{P}$ 800 MHz to 2.5 GHz</p> <p>Where P is the maximum output power rating of the transmitter in watts [W] according to the transmitter manufacturer and d is the recommended separation in meters [m].</p> <p>Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey^a, should be less than the compliance level in each frequency range^b.</p> <p>Interference may occur in the vicinity of equipment marked with the following symbol:</p> 
<p>Note: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection by structures, objects and people.</p>			
<p>^a Field strengths from fixed transmitters, such as base stations for radio [cellular/cordless] telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the unit is used exceeds the applicable RF compliance level above, the unit should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as reorienting or relocating the unit.</p> <p>^b Over the frequency range 150 kHz to 80 MHz, field strengths should be less than 3 V/m.</p>			

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Table 206 Recommended separation distances between portable and mobile RF communications equipment and the HALOGEN 250 twin model 201133 20			
The HALOGEN 250 twin model 201133 20 is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or user of the unit can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the unit as recommended below, according to the maximum output power of the communications equipment.			
Rated maximum output power of transmitter [W]	Separation distance d [m] according to frequency of transmitter		
	150 kHz to < 80 MHz $d = \left[\frac{3.5}{3} \right] \sqrt{P}$	80 MHz to < 800 MHz $d = \left[\frac{3.5}{3} \right] \sqrt{P}$	800 MHz to 2.5 GHz $d = \left[\frac{7}{3} \right] \sqrt{P}$
0.01	0.12	0.12	0.23
0.1	0.38	0.38	0.73
1	1.2	1.2	2.3
10	3.8	3.8	7.3
100	12	12	23
For transmitters rated at a maximum output power not listed above, the recommended separation distance d in meters (m) can be estimated using the equation applicable to the frequency of the transmitter, where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.			
Note: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.			

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